

**AVAILABLE IN:
CUSTOM SIZES • HEIGHTS • CONFIGURATIONS • DESIGNS**

WHY THE EURO?

If cosmetics are important in your mini-environment, consider the myriad advantages of the EURO Design Blower Filter Module. This unit offers a rivetless, low-profile design with a host of options to meet your critical needs.



FEATURES

- Rivetless Design
- Low Profile - 9 inch
- Backward Curved Motorized Impeller
(EBM Motor: CE • UL • CSA)
- Aluminum Construction
- Anodized Finish
- Sealant - Silicone Free
- Poron Dovetail Gasket
- Speed Control
- Integral Latches



OPTIONS

- Electrical:
 - EC
 - 230 VAC
 - 115 VAC
- Power Inlet:
 - IEC
 - Power Cord
 - NEMA Twist Lock
 - Amp
 - Harting
 - Conduit Whip

Filter Media:

- HEPA or ULPA
- Microfine Fiberglass: Pack Depth 2", 3" and 4"
- PTFE: 99.99995% @ MPPS for use in 300mm SEMI specifications

Filter Seals:

- Knife Edge
- Gel

Filter Screen:

- Expanded Painted Metal (*standard*)
- Expanded Anodized Aluminum
- Expanded Stainless Steel
- Perforated Stainless Steel

Also Available:

- Differential Pressure Sensors
- Power Indicator Light
- Motor Running Lights
- Shroud for Lighting & Ionization

DISTRIBUTED BY